

# ECG Noise and Lead-Failure Stress-Test Corpus (v1) — Generation Report and Full Record Roster

Ioannis Valasakis

2026-07-13

## ECG Noise and Lead-Failure Stress-Test Corpus (v1)

Ioannis Valasakis

*Electrocardiography Group, University of Glasgow*

**Date:** 2026-07-13 · **Repository:** <https://github.com/depolarised/artefaux> (tag v1.0.0) · **Sources** (every value below traces to these): `recipes/corpus.yaml`, `recipes/source_ids/*.csv`, `manifest.json` (all committed), regenerated in-memory from `artefaux.corpus + artefaux.recipes` (master\_seed = 20260713). Generator: `scripts/reports/generate_c`

**Scope.** This is a *stress-test set*, not a *clinically representative cohort*. It deliberately over-samples failure; do not use it to estimate real-world prevalence or deployment accuracy.

---

### 1. Objectives

This is a deterministic generator and reproducible corpus **definition** for stress-testing ECG signal-quality gates (the internal **signalguard**) and noise detectors (**noiseguard**). It does **not** redistribute derived signals: it ships code, recipes, labels, manifest, and provenance so a user regenerates the corpus **bit-exactly** from their own open PhysioNet copies.

This report documents the generation algorithm and enumerates **all 67 stress records** with their corruption truth and expected behaviour. Records resolve against a local PhysioNet mirror; the concrete PTB-XL parent ids shown here are those in `recipes/source_ids/` (reproducible, seeded selection).

---

### 2. Corpus composition

| Group          | n  | Source                                          | Corruption              |
|----------------|----|-------------------------------------------------|-------------------------|
| Naturally poor | 15 | PTB-XL<br>technical-validation<br>quality flags | none (inherently noisy) |

| Group                | n         | Source                                     | Corruption                               |
|----------------------|-----------|--------------------------------------------|------------------------------------------|
| Real-noise pairs     | 30        | clean PTB-XL + NSTDB <code>em/ma/bw</code> | SNR ladder $\{-6, +0, +6, +12, +18\}$ dB |
| Engineering extremes | 22        | clean PTB-XL parents (+ MACECGDB motion)   | single-lead (12) + multi-lead (10)       |
| <b>Total</b>         | <b>67</b> |                                            | 11 labelled data-integrity failures      |

Each record carries three label layers — **clinical parent** (authored rhythm class + PTB-XL quality flags; the Uni-G statement field is reserved and empty in v1), **corruption truth** (electrodes/leads, requested + measured SNR, noise segment, seed,  $\alpha$ , amplitude bookkeeping), and **expected behaviour** (per-lead and record-level, mapped to the `signalguard` `QualityLabel` / `RecordQuality` and `noiseguard` discard vocabularies).




Figure 1: Corpus composition — 67 stress records plus paired clean parents.

### 3. Sources and licensing

| Dataset       | Role                                                              | Version | License    |
|---------------|-------------------------------------------------------------------|---------|------------|
| PTB-XL        | clean + naturally-poor parents (500 Hz)                           | 1.0.3   | CC-BY-4.0  |
| PTB-XL+       | Uni-G/Glasgow layer — <i>cited and pinned, not consumed by v1</i> | 1.0.1   | CC-BY-4.0  |
| MIT-BIH NSTDB | <code>em/ma/bw</code> noise for the SNR ladder                    | 1.0.0   | ODC-BY-1.0 |
| MACECGDB      | real standing/walking/jumping motion                              | 1.0.0   | ODC-BY-1.0 |

PTB-XL, PTB-XL+, and MACECGDB are read from a local PhysioNet mirror; **only NSTDB** is fetched by `make download`. **Code** is GPL-3.0-or-later; **manifest, labels, and docs** are CC-BY-4.0. No source signals are redistributed — see `ATTRIBUTION.md`.

## 4. Generation algorithm

**4.1 Canonical form.** Every parent becomes a (12, T) float64 mV array in canonical lead order [I, II, III, aVR, aVL, aVF, V1-V6] at 500 Hz (NSTDB resampled from 360 Hz by deterministic polyphase resampling). Absent limb leads are derived from I and II.

**4.2 Real-noise SNR mixing.** For a clean lead  $x$ , real noise  $n$ , and requested SNR  $s$  (dB):  $\alpha = \sqrt{P_x / (P_n \cdot 10^{s/10})}$ ,  $y = x + \alpha n$ , with  $P_x$ ,  $P_n$  the DC-removed mean-square powers over the contaminated interval only. Both the requested and the *measured* SNR, and the exact NSTDB segment offset, are recorded in the label.




Figure 2: Real-noise SNR ladder (illustrative synthetic parent; NSTDB-style electrode motion).

**4.3 Electrode-domain corruption.** Limb-electrode potentials are recovered in the RA-as-reference gauge (RA = 0, LA = I, LL = II); chest-electrode potentials as  $Ck = Vk + WCT$ ,  $WCT = (RA+LA+LL)/3$ . An artefact is added to a chosen electrode; the WCT and all twelve leads are recomputed, and the invalidated derived leads are recorded. (The classical WCT is an idealised electrode model, **not** a torso forward model.)

**4.4 Engineering extremes.** Deterministic builders for clipping/rail-saturation, flatline (zero), constant-ADC, lead-off, digital-missing (NaN), electrode-motion swings, step-and-exponential-

recovery, polarity inversion, and intermittent lead-off. NaN / flat / constant / rail cases are labelled `data_integrity_failure` with a type — **not** as noise.



Figure 3: Engineering extreme: an overload swing clipped at the acquisition rail.

## 5. Full record roster

All 67 records are listed below. **Expected quality** is the record-level `RecordQuality` a correctly-tuned gate should assign; **Discard** is the expected `noiseguard` record-level decision; **Bad leads** is the count of leads the label expects a detector to reject. PTB-XL parents are shown as record stems (the `records500/NNNNN/` directory prefix is omitted).

### 5.1 Naturally poor (15)

Real PTB-XL records selected by technical-validation quality flags; no clean parent (inherently bad).

| Record                   | PTB-XL source         | Quality flag(s)                  | Expected quality | Discard |
|--------------------------|-----------------------|----------------------------------|------------------|---------|
| <code>nlf_nat_001</code> | <code>10646_hr</code> | <code>electrodes_problems</code> | <b>reject</b>    | yes     |
| <code>nlf_nat_002</code> | <code>10652_hr</code> | <code>burst_noise</code>         | <b>limited</b>   | —       |
| <code>nlf_nat_003</code> | <code>11176_hr</code> | <code>static_noise</code>        | <b>limited</b>   | —       |
| <code>nlf_nat_004</code> | <code>12587_hr</code> | <code>static_noise</code>        | <b>limited</b>   | —       |
| <code>nlf_nat_005</code> | <code>13777_hr</code> | <code>static_noise</code>        | <b>limited</b>   | —       |
| <code>nlf_nat_006</code> | <code>14197_hr</code> | <code>static_noise</code>        | <b>limited</b>   | —       |
| <code>nlf_nat_007</code> | <code>16149_hr</code> | <code>baseline_drift</code>      | <b>reject</b>    | yes     |
| <code>nlf_nat_008</code> | <code>16981_hr</code> | <code>static_noise</code>        | <b>limited</b>   | —       |
| <code>nlf_nat_009</code> | <code>17281_hr</code> | <code>static_noise</code>        | <b>limited</b>   | —       |
| <code>nlf_nat_010</code> | <code>17293_hr</code> | <code>baseline_drift</code>      | <b>limited</b>   | —       |
| <code>nlf_nat_011</code> | <code>19487_hr</code> | <code>static_noise</code>        | <b>limited</b>   | —       |
| <code>nlf_nat_012</code> | <code>20478_hr</code> | <code>static_noise</code>        | <b>limited</b>   | —       |
| <code>nlf_nat_013</code> | <code>00359_hr</code> | <code>static_noise</code>        | <b>reject</b>    | yes     |
| <code>nlf_nat_014</code> | <code>06320_hr</code> | <code>baseline_drift</code>      | <b>limited</b>   | —       |
| <code>nlf_nat_015</code> | <code>06968_hr</code> | <code>static_noise</code>        | <b>limited</b>   | —       |

## 5.2 Real-noise SNR-ladder pairs (30)

Clean PTB-XL parents mixed with NSTDB noise across the ladder {-6, +0, +6, +12, +18} dB. Each ships with its untouched clean parent (\*\_clean).

| Record        | Rhythm          | PTB-XL parent | Noise @ SNR | Leads | Bad leads | Expected quality  | Discard |
|---------------|-----------------|---------------|-------------|-------|-----------|-------------------|---------|
| nlf_noise_001 | sinus           | 10236_hr      | em @ -6 dB  | 12    | 12        | <b>reject</b>     | yes     |
| nlf_noise_002 | afib            | 10330_hr      | ma @ +0 dB  | 6     | 0         | <b>limited</b>    | yes     |
| nlf_noise_003 | pac pvc         | 11089_hr      | bw @ +6 dB  | 6     | 0         | <b>limited</b>    | —       |
| nlf_noise_004 | bbb con-duction | 11662_hr      | em @ +12 dB | 12    | 0         | <b>diagnostic</b> | —       |
| nlf_noise_005 | st t            | 12086_hr      | ma @ +18 dB | 6     | 0         | <b>diagnostic</b> | —       |
| nlf_noise_006 | flutter svt     | 12407_hr      | bw @ -6 dB  | 6     | 6         | <b>limited</b>    | yes     |
| nlf_noise_007 | sinus           | 01306_hr      | em @ +0 dB  | 12    | 0         | <b>limited</b>    | yes     |
| nlf_noise_008 | afib            | 13114_hr      | ma @ +6 dB  | 6     | 0         | <b>limited</b>    | —       |
| nlf_noise_009 | pac pvc         | 13332_hr      | bw @ +12 dB | 6     | 0         | <b>diagnostic</b> | —       |
| nlf_noise_010 | bbb con-duction | 13611_hr      | em @ +18 dB | 12    | 0         | <b>diagnostic</b> | —       |
| nlf_noise_011 | st t            | 14071_hr      | ma @ -6 dB  | 6     | 6         | <b>limited</b>    | yes     |
| nlf_noise_012 | flutter svt     | 14211_hr      | bw @ +0 dB  | 6     | 0         | <b>limited</b>    | yes     |
| nlf_noise_013 | sinus           | 15473_hr      | em @ +6 dB  | 12    | 0         | <b>limited</b>    | —       |
| nlf_noise_014 | afib            | 15535_hr      | ma @ +12 dB | 6     | 0         | <b>diagnostic</b> | —       |
| nlf_noise_015 | pac pvc         | 15877_hr      | bw @ +18 dB | 6     | 0         | <b>diagnostic</b> | —       |

| Record                | Rhythm              | PTB-XL<br>parent | Noise @<br>SNR | Leads | Bad leads | Expected<br>quality | Discard |
|-----------------------|---------------------|------------------|----------------|-------|-----------|---------------------|---------|
| nlf_<br>noise_<br>016 | bbb con-<br>duction | 16615_hr         | em @ -6<br>dB  | 12    | 12        | <b>reject</b>       | yes     |
| nlf_<br>noise_<br>017 | st t                | 16930_hr         | ma @ +0<br>dB  | 6     | 0         | <b>limited</b>      | yes     |
| nlf_<br>noise_<br>018 | flutter<br>svt      | 17553_hr         | bw @ +6<br>dB  | 6     | 0         | <b>limited</b>      | —       |
| nlf_<br>noise_<br>019 | sinus               | 17727_hr         | em @ +12<br>dB | 12    | 0         | <b>diagnostic</b>   | —       |
| nlf_<br>noise_<br>020 | afib                | 17770_hr         | ma @ +18<br>dB | 6     | 0         | <b>diagnostic</b>   | —       |
| nlf_<br>noise_<br>021 | pac pvc             | 17867_hr         | bw @ -6<br>dB  | 6     | 6         | <b>limited</b>      | yes     |
| nlf_<br>noise_<br>022 | bbb con-<br>duction | 18284_hr         | em @ +0<br>dB  | 12    | 0         | <b>limited</b>      | yes     |
| nlf_<br>noise_<br>023 | st t                | 18508_hr         | ma @ +6<br>dB  | 6     | 0         | <b>limited</b>      | —       |
| nlf_<br>noise_<br>024 | flutter<br>svt      | 19002_hr         | bw @ +12<br>dB | 6     | 0         | <b>diagnostic</b>   | —       |
| nlf_<br>noise_<br>025 | sinus               | 19238_hr         | em @ +18<br>dB | 12    | 0         | <b>diagnostic</b>   | —       |
| nlf_<br>noise_<br>026 | afib                | 01979_hr         | ma @ -6<br>dB  | 6     | 6         | <b>limited</b>      | yes     |
| nlf_<br>noise_<br>027 | pac pvc             | 20943_hr         | bw @ +0<br>dB  | 6     | 0         | <b>limited</b>      | yes     |
| nlf_<br>noise_<br>028 | bbb con-<br>duction | 20977_hr         | em @ +6<br>dB  | 12    | 0         | <b>limited</b>      | —       |
| nlf_<br>noise_<br>029 | st t                | 21495_hr         | ma @ +12<br>dB | 6     | 0         | <b>diagnostic</b>   | —       |
| nlf_<br>noise_<br>030 | flutter<br>svt      | 21820_hr         | bw @ +18<br>dB | 6     | 0         | <b>diagnostic</b>   | —       |

### 5.3 Engineering extremes — single-lead (12)

| Record                              | PTB-XL<br>parent | Corruption                                      | Integrity<br>failure    | Bad leads | Expected<br>quality | Discard |
|-------------------------------------|------------------|-------------------------------------------------|-------------------------|-----------|---------------------|---------|
| nlf_eng_001_emg_v1                  | 02965_hr         | swing 6 mV p2p (V1)                             | —                       | 1         | <b>limited</b>      | —       |
| nlf_eng_002_macecg_walk_v2          | 03016_hr         | MACECGDB walk swing 8 mV p2p (V2)               | —                       | 1         | <b>limited</b>      | —       |
| nlf_eng_003_intermittent_v6         | 03047_hr         | intermittent lead-off with reconnection (V6)    | intermittent lead off   | 1         | <b>limited</b>      | —       |
| nlf_eng_004_step_ii                 | 03365_hr         | baseline step 4 mV, tau=1 s (II)                | —                       | 0         | <b>limited</b>      | —       |
| nlf_eng_005_macecg_jump_overload_v3 | 03427_hr         | MACECGDB jump swing 30 mV p2p (V3) → rail ±5 mV | rail saturation         | 1         | <b>limited</b>      | —       |
| nlf_eng_006_flatline_v4             | 03543_hr         | flatline / zero-fill (V4)                       | flatline zero           | 1         | <b>limited</b>      | —       |
| nlf_eng_007_constant_v5             | 03687_hr         | stuck constant 1.2 mV (V5)                      | constant adc            | 1         | <b>limited</b>      | —       |
| nlf_eng_008_missing_v6              | 03730_hr         | digital-missing / NaN (V6)                      | digital missing channel | 1         | <b>limited</b>      | —       |
| nlf_eng_009_flatline_i              | 03941_hr         | flatline / zero-fill (I)                        | flatline zero           | 5         | <b>limited</b>      | —       |
| nlf_eng_010_leadoff_v2              | 04593_hr         | lead-off (V2)                                   | lead off                | 1         | <b>limited</b>      | —       |
| nlf_eng_011_emg_burst_v5            | 00476_hr         | swing 5 mV p2p (V5), t in [3,6] s               | —                       | 0         | <b>limited</b>      | —       |
| nlf_eng_021_reversal_i              | 09509_hr         | polarity inversion (I)                          | —                       | 0         | <b>limited</b>      | —       |

## 5.4 Engineering extremes — multi-lead / mixed (10)

Electrode-domain propagation, precordial coupling, multi-lead saturation, and compound “wild” cases.

| Record                                  | PTB-XL<br>parent | Corruption                                                                                                | Integrity<br>failure | Bad leads | Expected<br>quality | Discard |
|-----------------------------------------|------------------|-----------------------------------------------------------------------------------------------------------|----------------------|-----------|---------------------|---------|
| nlf_eng_<br>012_<br>couple_<br>v1v2     | 05267_hr         | precordial<br>coupling<br>V1, V2 @<br>+0 dB                                                               | —                    | 2         | <b>limited</b>      | —       |
| nlf_eng_<br>013_<br>couple_<br>v4v5     | 05298_hr         | precordial<br>coupling<br>V4, V5 @<br>+6 dB                                                               | —                    | 0         | <b>limited</b>      | —       |
| nlf_eng_<br>014_<br>triple_<br>v1v2v3   | 05887_hr         | precordial<br>coupling<br>V1, V2, V3<br>@ -6 dB                                                           | —                    | 3         | <b>limited</b>      | yes     |
| nlf_eng_<br>015_<br>electrode_<br>1a    | 06175_hr         | LA<br>electrode<br>motion 1<br>mV                                                                         | —                    | 4         | <b>limited</b>      | yes     |
| nlf_eng_<br>016_<br>electrode_<br>ra    | 06526_hr         | RA<br>electrode<br>offset 3 mV                                                                            | —                    | 6         | <b>reject</b>       | yes     |
| nlf_eng_<br>017_<br>saturate_<br>v2v3   | 06789_hr         | overload<br>swing 30<br>mV p2p<br>(V2, V3) →<br>rail ±5 mV                                                | rail<br>saturation   | 2         | <b>limited</b>      | —       |
| nlf_eng_<br>018_<br>saturate_<br>v4v5v6 | 07086_hr         | overload<br>swing 30<br>mV p2p<br>(V4, V5,<br>V6) → rail<br>±5 mV                                         | rail<br>saturation   | 3         | <b>limited</b>      | yes     |
| nlf_eng_<br>019_<br>electrode_<br>1l    | 08699_hr         | LL<br>electrode<br>motion 1.2<br>mV                                                                       | —                    | 3         | <b>limited</b>      | yes     |
| nlf_eng_<br>020_<br>compound_<br>wild   | 08931_hr         | LA<br>electrode<br>motion 0.8<br>mV ;<br>MACECGDB<br>walk swing<br>8 mV p2p<br>(V2) ;<br>lead-off<br>(V5) | lead off             | 5         | <b>limited</b>      | yes     |



| Record                                  | PTB-XL<br>parent | Corruption                                                                                                                                     | Integrity<br>failure          | Bad leads | Expected<br>quality | Discard |
|-----------------------------------------|------------------|------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------|-----------|---------------------|---------|
| nlf_eng_<br>022_<br>mixed_<br>integrity | 09915_hr         | digital-<br>missing /<br>NaN (V1) ;<br>overload<br>swing 25<br>mV p2p<br>(V6) → rail<br>±5 mV ;<br>baseline<br>step 3 mV,<br>tau=1.5 s<br>(II) | digital<br>missing<br>channel | 2         | <b>limited</b>      | —       |

## 6. Reproducibility and provenance

A single master seed (20260713) plus `provenance.json` (pinned source versions, per-record seed, source-file SHA-256) make the corpus bit-exactly regenerable:

```
make download      # NSTDB only (PTB-XL / PTB-XL+ / MACECGDB from the local mirror)
make select        # resolve PTB-XL parent ids -> recipes/source_ids/ (seeded)
make regenerate    # WFDB records + labels + manifest + provenance -> ./out
```

Determinism is per-record: `SeedSequence(entropy=master_seed, spawn_key=(record_index,))`. The corpus definition leaves source ids unresolved; `resolve_source_ids()` binds the seeded PTB-XL selection into the specs at generation time.

## 7. Verification

Confirmed end-to-end on the local PhysioNet mirror (PTB-XL 1.0.3, MACECGDB 1.0.0, NSTDB 1.0.0):

- **Counts:** 67 records (15 / 30 / 22) + 52 paired clean parents.
- **SNR ladder:** max |requested - measured| SNR error **0.0 dB** across all 30 mixes (tol 0.5 dB).
- **Data integrity:** digital-missing (NaN) round-trips through WFDB — the affected channel reads back all-NaN while intact leads stay finite.
- **Determinism:** two independent full regenerations are **byte-identical**.
- **Quality gates:** 61 unit/integration tests pass; `ruff` + `black` clean.

Generated by Artefaux v1.0.0 (<https://github.com/depolarised/artefaux>) via `scripts/reports/generate_corpus_` every figure and table value is computed from the committed corpus definition.